

June 2015

Subject: Broadley-James D.O. Membrane Cartridge
Model Numbers KA*, KB*, KC*, MS*, MK*

To Whom It May Concern:

The product contact polymer is Silicone; a backing layer of Teflon is not a wetted material and does not contact product. The polymers used in the membrane are manufactured in accordance with applicable portions of FDA CFR 117.1550 and CFR 177.2600. The Silicone formula has been tested and passed USP Class VI testing at ~70°C.

Broadley-James Corporation is dedicated to quality and maintains a high degree of control over all aspects of production. The products we ship are individually pressure tested, measured, and compared to Broadley-James' internal manufacturing standards, and must meet or exceed those standards.

Thank you for choosing Broadley-James Corporation.

Best Regards,



Robert Garrahy
Vice President Bioprocess Technologies

/ses