

How do I choose the right membrane cartridge for my sensor?

Membrane Cartridge Selection		
If your sensor cathode looks like this		
<p>Broadley James OxyProbe II</p> 	<p>Broadley James OxyProbe</p> 	<p>Broadley James OxyProbe</p> 
<p>Mettler Toledo® Ingold®</p> 	<p>Mettler Toledo® Ingold®</p> 	<p>Mettler Toledo® Ingold®</p> 
Use this membrane cartridge model		
<p>OxyProbe II</p> 	<p>OxyProbe (12mm)</p> 	<p>OxyProbe (25mm)</p> 

The contents of this publication are presented for information purposes only, and while every effort has been made to ensure their accuracy, they are not to be construed as warranties or guarantees, express or implied, regarding products or services described herein or in their use or applicability. All sales are governed by our terms and conditions, which are available on request. We reserve the right to modify or improve the design or specification of such products at any time.

© 2021 Broadley-James Corporation. All rights reserved. Visit www.broadley-james.com/trademarks for trademark information.

TMP-BF-102101